

# Sani Wipes Lemon



Hand sanitizer wipes



**SANI WIPES LEMON** hand sanitizer wipes are a one-step system combining both an abrasive non-scratching wipe and a highly effective skin sanitizing formula. **SANI WIPES LEMON** hand sanitizer wipes absorb and remove disease causing germs and bacteria. This product contains a highly effective antimicrobial formula that assures a rapid kill. **SANI WIPES LEMON** hand sanitizer wipe's unique formula reduces the risk of infection by killing 99.9% of disease causing organisms.

Micro-Organism (Exposure Time):

- Pseudomonas aeruginosa (10 Seconds) • Staphylococcus aureus (10 Seconds) • Escherichia Coli (20 Seconds) • Streptococcus pyogenes (10 Seconds) • Salmonella choleraesuis (10 Seconds) • Candida albicans (10 Seconds) • Staphylococcus aureus - MRSA (10 Seconds).

**Directions:** Use to supplement regular handwashing. Remove lid. Pull towels up from the center of the roll and twist into a point. Push twisted towel through center of notched opening in the underside of the lid. Replace lid, making sure lid fits securely. Pull out wipe at a 90 degree angle and tear off at perforation. Dispose of properly. Remaining wipes feed automatically for next use. Always reseal cap to insure against moisture loss and contamination.

Health	1
Flammability	2
Reactivity	0
Protection	A

#### Areas of Use:

- \*school districts
- \*cities and counties
- \*hospitals
- \*offices
- \*janitorial services
- \*retirement facilities
- \*hotels

Appearance and odor.....	Saturated towel, Lemon Scent
Specific Gravity.....	Concentrate - 0.97
VOC.....	N/A
NFPA Fire Rating.....	Flammable

## ONLY OMEGA™

Fairfield, CA 94534

1-800-571-7347 [www.onlyomega.com](http://www.onlyomega.com) Fax: 707-864-8134